

TNamed
fName fTitle fglsA
@~TNamed TNamed TNamed TNamed operator= Clear Clone Compare Copy FillBuffer GetName GetTitle Hash IsSortable SetName SetNameTitle SetTitle Is Print Sizeof Class Class_Name IsA ShowMembers

HParSet
detName[10] versions[3] status changed paramContext author description fglsA
@~HParSet HParSet init init init write write clear print getDetectorName resetInputVersions setInputVersion getInputVersion setStatic isStatic setChanged hasChanged setParamContext getParamContext setAuthor getAuthor setDescription getDescription copyComment Class Class_Name IsA ShowMembers Streamer StreamerNVirtual

HParCond
fglsA
@~HParCond putParams getParams printParams Class Class_Name IsA ShowMembers Streamer

HKickTrackFPar
fPMin fPMin2 fMaxXPull1 fMaxXPull2 fXPullStep fTofXResolution fMinTof fMaxTof fglsA
@~HKickTrackFPar HKickTrackFPar clear init write writeFile read getPMin getPMin2 getMaxXPull1 getMaxXPull2 getXPullStep getTofXResolution getMinTof getMaxTof printTo readFrom readFrom putParams getParams Class Class_Name IsA ShowMembers Streamer StreamerNVirtual