

TNamed
fName fTitle fglsA
@~TNamed TNamed TNamed TNamed TNamed operator= Clear Clone Compare Copy FillBuffer GetName GetTitle Hash IsSortable SetName SetNameTitle SetTitle Is Print Sizeof Class Class_Name IsA ShowMembers

HParSet
detName[10] versions[3] status changed paramContext author description fglsA
@~HParSet HParSet init init init write write clear print getDetectorName resetInputVersions setInputVersion getInputVersion setStatic isStatic setChanged hasChanged setParamContext getParamContext setAuthor getAuthor setDescription getDescription copyComment Class Class_Name IsA ShowMembers Streamer StreamerNVirtual

HRichParSet
fglsA m_pReadParam
@~HRichParSet HRichParSet init write Class Class_Name IsA ShowMembers Streamer StreamerNVirtual

HRichMappingPar
uiplug padsxy padsx padsy xyuip xyuiplog uncuip i fPadsNr fCols fRCs fPorts fModules fChannels fRows fAddress1 fAddress2 fUnconnCh fglsA
@~HRichMappingPar HRichMappingPar initAscii writeAscii getRC getPort getMod getCh getCh isReadOut getCol getRow isValidUpi isUnConnCh Class Class_Name IsA ShowMembers Streamer StreamerNVirtual