

TNamed
fName fTitle fgIsA
@~TNamed TNamed TNamed TNamed TNamed operator= Clear Clone Compare Copy FillBuffer GetName GetTitle Hash IsSortable SetName SetNameTitle SetTitle Is Print Sizeof Class Class_Name IsA ShowMembers

HParSet
detName[10] versions[3] status changed paramContext author description fgIsA
@~HParSet HParSet init init init write write clear print getDetectorName resetInputVersions setInputVersion getInputVersion setStatic isStatic setChanged hasChanged setParamContext getParamContext setAuthor getAuthor setDescription getDescription copyComment Class Class_Name IsA ShowMembers Streamer StreamerNVirtual

HRichParSet
fgIsA m_pReadParam
@~HRichParSet HRichParSet init write Class Class_Name IsA ShowMembers Streamer StreamerNVirtual

HRichSectorPar
fRotationAngle fSectorId fgIsA
@~HRichSectorPar HRichSectorPar reset setSectorId getSectorId getPhi setPhi Class Class_Name IsA ShowMembers Streamer StreamerNVirtual